

IHE Integration Statement			
Vendor	Product Name	Version	Date
Quadramed	Smart Identity Exchange	2012.1	January 13, 2012
This product implement all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
Profiles Implemented	Actors Implemented	Options Implemented	
Audit Trail and Node Authentication	Secure Application	None	
Consistent Time	Time Client		
Patient Demographics Query	Patient Demographics Supplier	Pediatric Demographics	
Patient Demographics Query HL7 V3	Patient Demographics Supplier	Continuation Option, Pediatric Demographics	
Patient Identifier Cross-referencing for MPI	Patient Identity Cross-reference Manager	Pediatric Demographics	
Patient Identifier Cross-Reference HL7 V3	Patient Identity Cross-reference Manager	Pediatric Demographics	
<b>Internet address for vendor's IHE information:</b>			
<a href="http://www.quadramed.com/IHE2012.pdf">http://www.quadramed.com/IHE2012.pdf</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			
<b>HL7:</b> <a href="http://hl7.org/">http://hl7.org/</a>			
<b>DICOM:</b> <a href="http://medical.nema.org/">http://medical.nema.org/</a>			
<b>Links to general information on IHE</b>			
North America: <a href="http://www.rsna.org/IHE">http://www.rsna.org/IHE</a>	Europe: <a href="http://www.ihe-europe.org">http://www.ihe-europe.org</a>	Japan: <a href="http://www.ihe-j.org">www.ihe-j.org</a>	